


QuickPrep™ Tips

Fast and easy to use extraction pipette tip

QuickPrep tips save money and increase lab throughput

QuickPrep tips, powered by ABS Sorbent technology, are an ideal yes/no indicator of the presence of compounds and can be used for the quick purification and separation of samples prior to analysis. QuickPrep tips are fast and easy to use.



Benefits of QuickPrep™ Tips

- **Fast sample clean-up** for LC-MS and GC-MS
- Concentrate analytes in **less than 1 minute**
- Use in manual or automated pipetting systems
- Extract a **wide range of analytes** - from polar to non-polar compounds
- Faster and less expensive than solid phase extraction (SPE)

For a list of compounds captured by QuickPrep tips, please visit: www.absmaterials.com/technical-resources.aspx

Fast

QuickPrep Tips extract a wide range of analytes and concentrate the analytes from aqueous samples into organic solvents in less than one minute.

Easy

Processing samples by pipetting with QuickPrep Tips is much easier and less time consuming than running traditional SPE methods.

Affordable

Less expensive than traditional SPE methods, QuickPrep tips are priced affordably, saving money and increasing lab throughput.

QuickPrep Tips have a sample loading volume of **1000 uL (1mL)** and are available in **packs of 96**.

Purchase on Amazon.com, or on our webstore at: <http://www.absmaterials.com/abs-store.aspx>